
Phosphate Bonded Investments


Ceramigold

Investment with carbon for ceramic gold and crown and bridge alloys

Revestimiento con carbón para aleaciones de oro cerámico, coronas y puentes
Rivestimento contenente grafite per leghe oro ceramica e leghe per ponti e corone
Revêtement graphité pour or céramique et alliages pour couronne et bridge
Kohlenstoffhaltige Einbettmasse für Goldkeramik sowie Kronen- und Brückenlegierungen

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
Hi-Temp

Carbon-free investment for high fusing alloys

Revestimiento sin carbón para aleaciones de alto grado de fusión
Rivestimento non contenente grafite per leghe ad alto punto di fusione
Revêtement sans graphite pour des alliages à point de fusion élevé
Kohlenstofffreie Einbettmasse für hochschmelzende Legierungen

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Cera-Fina

Carbon-free fine grained investment for high fusing alloys

Revestimiento de grano fino sin carbón para aleaciones de alto grado de fusión
Rivestimento a grana fine, non contenente grafite, per leghe ad alto punto di fusione
Revêtement à granulation fine sans graphite pour des alliages à point de fusion élevé

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1. Pattern Preparation:

- Attach sprue to the wax pattern using pliable sticky wax.
- Mount patterns on crucible former base.
- Spray with Whip Mix SMOOTHX Wax Pattern Cleaner. Gently blow excess cleaner from pattern.

2. Ring Preparation:

- Line casting ring with liner positioned 6 mm short on each end.
- Tack down in place with pliable sticky wax.
- Soak lined ring in water for a minimum of one minute, then shake well (5 times).
- Place ring onto former base.

3. Mixing:

Suggested Liquid/Powder Ratio:

	60 g	90 g	100 g
Ceramigold	9.5 mL	14.5 mL	16 mL
Hi-Temp	9.5 mL	14.5 mL	16 mL
Cera-Fina	14.5 mL	21.5 mL	24 mL

- Rinse bowl out with water — shake out excess.
- Add measured liquid into bowl then add weighed powder.
- Hand mix well for 20 seconds.
- Mechanical mix under vacuum, slow speed (350–450 RPM) for 60–90 seconds.
- Benchset for a minimum of 1 hour. For large, ringless molds, bench set for a minimum of 2 hours. (**Cera-Fina**: Bench set for a minimum of 1½ hours. For large, ringless molds, bench set for a minimum of 2 hours.)

4. Burnout:

- Remove former base and trim glaze off top of mold.

When using all wax patterns and sprues:

- Place mold in furnace at room temperature and raise temperature at a rate of 8°C (15°F) per minute to 430°C (800°F).

- Hold for 30 minutes.
- Raise to maximum temperature at a rate of 14°C (25°F) per minute.
- Hold for 30 minutes.

If using plastic, large wax bridges, large ringless molds or precision attachments:

- Place mold in furnace at room temperature and raise temperature at a rate of 4°C (8°F) per minute to 430°C (800°F).
- Hold for 30 to 60 minutes.
- Raise to maximum temperature at a rate of 7°C (12°F) per minute.
- Hold for 30 minutes.

5. Casting:

- Cast alloy with oxygen/propane torch or induction machine according to alloy manufacturer's instructions.
- Clean casting with aluminum oxide (50–60 microns), fine blasting sand or glass beads.

Note:

- Recommended Liquid: Special Liquid Concentrate
- Use separate mixing bowls for phosphate and gypsum investments.
- For best results, keep the liquid and powder at room temperature of about 23°C (73°F).
- For ringless molds larger than 270 g, do not use liquid concentrations over 75%.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

Recommended Liquid Concentrations:

Alloy	Liquid Concentration	Burnout Temperature
Chrome-Cobalt	100%	870–900°C (1600–1650°F)
Nickel-Beryllium	75–100%	815–870°C (1500–1600°F)
Palladium (2–8% Gold)	10–15%	700–870°C (1300–1600°F)
Silver-Palladium	75%	760–815°C (1400–1500°F)
Ceramic Gold (45% or higher)	50–75%	700–815°C (1300–1500°F)
Crown & Bridge Gold	50%	575–700°C (1100–1300°F)

Physical Properties:

	Ceramigold	Hi-Temp	Cera-Fina
Liquid*/Powder Ratio	16 mL/100 g	16 mL/100 g	24 mL/100 g
Working Time	6–7 minutes	7–8 minutes	8–10 minutes
Setting Expansion	0.7%	0.7%	0.25%
Thermal Expansion	1.2%	1.2%	1.3%
Compressive Strength, wet	1,500 psi (10 MPa)	1,500 psi (10 MPa)	600 psi (4 MPa)

*Suggested concentration of Special Liquid is 75% (3 parts liquid to 1 part water).

** Recommended concentrations are approximate and can be adjusted to optimize fit. To increase expansion, use more liquid and less water. To decrease expansion, use less liquid and more water. Always maintain total liquid / water volume. For expansion higher than results achieved with Special Liquid Concentrate, use Special Liquid Concentrate-Plus

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