VERIGUIDE™ OS

3D PRINT RESIN
INSTRUCTIONS

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INDICATIONS FOR USE

Whip Mix VeriGuide™ OS is a 3D printed resin for producing drilling guides used in dental surgery. VeriGuide™ OS is formulated for 3D-printing systems (DLP, SLA, and LCD) at 385 nm and 405 nm. VeriGuide™ OS was successfully tested for biocompatibility and meets all mechanical properties requirement.

STORAGE

- VeriGuide™ OS must be stored in the original packaging at room temperature in a dry area.
- Always keep container tightly sealed and close the container immediately after each use.
- Do not expose VeriGuide™ OS to any light source. A minimal amount of light can induce polymerization (solidification).
- Resin that is stored in a tray must be stirred before each print using a plastic or paper card (similar to a business card), to ensure a homogenous mixture.
- Do not use the product past the expiration date.
- The lot number and expiration date are indicated on VeriGuide™ OS packaging. In case of questions or concerns, please refer to the lot number.

SURGICAL GUIDE DESIGN RECOMMENDATION

- Wall Thickness 2 mm
- Internal offset ≥ 0.1 mm
- Block out angle 0°
- Retention amount 0 mm
- Offset from Sleeve 0.10 mm

PROCESSING

- When printing with VeriGuide™ OS 3D Print Resin, make sure resin is at a temperature between 20°C and 30°C. Extreme high or low temperatures could potentially affect the accuracy of the printed models or cause failure of the printed objects.
- Please use VeriGuide™ OS material file in Asiga Composer software provided by Whip Mix Corporation.
- VeriGuide™ OS bottle must be shaken thoroughly for approximately five minutes before use.
- Surgical guide can be printed at slice thickness of 50, 75, or 100 microns.
- **Don't add supports to the intaglio surface.**
- Don't print with the end of the surgical guide arch attached to the Asiga build plate with any print orientation because failure could happen at the tip of the arch.
- After the build process is finished, post processing is required.

POST PROCESSING

- When the print job is complete, remove the build platform from the printer, then carefully remove the printed objects from the platform.
- **The printed object must be rinsed using a fresh ultrasonic alcohol bath for 5 – 10 minutes to eliminate any excess resin.**
- **Printed objects should then be rinsed in a second fresh ultrasonic alcohol bath for an additional 5 – 10 minutes.**
- To remove excess alcohol, leave prints to dry at room temperature (or compressed air can also be used).
- Carefully remove all support structures from the printed object.


- **Place the printed objects in the recommended light curing unit, Otoflash G171, for 6000 flashes or any equivalent curing unit with wavelength between 300 – 700 nm, at 200 Watt output, and for 10 minutes.**
- Look for cured debris in the tray after each print job. A paint strainer can be used to filter out any debris. Printing with cured debris could potentially cause damage to the tray or the printed objects.
- Use crosscut carbide bur at a slow speed to remove the remaining printed support pins
- Surgical guide can be disinfected and is compatible with sterilization methods such as steam/autoclave (up to 135°C for 5 minutes), ethylene oxide, and gamma irradiation. Be sure to follow sterilization guidelines. Please note that gamma irradiation may cause discoloration.

PHYSICAL PROPERTIES

TECHNICAL DATA:

- **Color:**
Clear
- **Density:**
ca. 1.1g/ml
- **Viscosity:** Average:
600 – 900 cP at 25°C

- **Storage:**

ORDERING INFORMATION:

Standard packing: 1kg (910ml)
Item no.: **71175**

| PROPERTY | RESULT | STANDARD |
|----------------------------|------------|-------------|
| Ultimate Flexural Strength | > 100 MPa | ISO 20795-2 |
| Flexural Modulus | > 2000 MPa | ISO 20795-2 |
| Residual Monomer | < 0.10 % | ISO 20795-2 |

Warning!

Combustible liquid
May cause an allergic skin reaction.
Suspected of damaging fertility.
May cause long lasting harmful effects to aquatic life.

Precautionary statements:

Obtain special instructions before use.
Do not handle until all safety precautions have been read and understood.
Keep away from flames and hot surfaces. No smoking.
Avoid breathing mists, vapors or spray.
Wash thoroughly after handling.

Contaminated work clothing must not be allowed out of the workplace.
Avoid release to the environment.
Wear protective gloves.
IF ON SKIN: Wash with plenty of soap and water.
If skin irritation or rash occurs: Get medical attention.
Take off contaminated clothing and wash it before reuse.
IF exposed or concerned: Get medical attention.
In case of fire: Use dry chemical, carbon dioxide or water spray.
Store in a well-ventilated place. Keep cool.
Store locked up.
Dispose of contents and container in accordance with local and national regulations.

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Whip Mix Corporation ■ 361 Farmington Avenue ■ Louisville, KY USA 40209