VERIMODEL™ OS

3D PRINT RESIN
INSTRUCTIONS

3D-DRUCKHARZ
ANLEITUNG

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PRODUCT DESCRIPTION

Whip Mix VeriModel™ OS is a 3D print resin for producing dental models on DLP, LCD or SLA technology printers with 385nm and 405nm light sources.

STORAGE

- VeriModel™ OS must be stored in the original packaging at room temperature in a dry area.
- Always keep container tightly sealed and close the container immediately after each use.
- Do not expose VeriModel™ OS to any light source. Minimum amount of light can induce polymerization (solidification).
- Once resin is poured in a tray, it should be stored in a dark area
- **Since pigments settle over time, resin that is stored in a tray must be stirred before each print using a plastic or paper card (similar to a business card) to ensure a homogenous mixture and to get pigments back into the solution. Failure to do this will obstruct the light, causing irregular edges / surfaces.**
- **NOTE: Avoid using sharp objects and edges as these can damage the resin tray.**
- Do not use the product past the expiration date.
- The lot number and expiration date are indicated on VeriModel™ OS packaging. In case of questions or concerns, please refer to the lot number.

PROCESSING

- When printing with VeriModel™ OS 3D Print Resin, make sure it is at a temperature between 20°C and 30°C. Extreme high or low temperatures could negatively affect the performance of the resin.
- VeriModel™ OS bottle must be shaken thoroughly for approximately five minutes before use. Color deviation may occur if not thoroughly shaken.
- Carefully pour VeriModel™ OS into the tray/vat.
- For optimized print settings refer to the Resource library at <https://whipmix.com/technical-resources/>
- After the build process is finished, post processing is required. If this cannot be immediately completed, leave the printed objects in the printer until ready to post process.

POST PROCESSING

- When the print job is complete, remove the build platform from the printer, then carefully remove the printed objects from the platform.
- Carefully remove all support structures from the printed object.
- **The printed objects must be rinsed using an ultrasonic alcohol bath for 3 – 5 minutes to eliminate any excess resin.**
- **Printed objects should then be rinsed in a second clean ultrasonic alcohol bath for an additional 5 minutes.**
- Models with cavities and gaps will take longer to clean.
- To remove excess alcohol, leave models to dry at room temperature (or compressed air can also be used).

- To achieve the maximum mechanical strength, post cure in a light curing unit, such as Asiga's Flash from Whip Mix. Follow manufacturer's instructions.
- **Look for cured debris in the tray after each print job, especially after a print failure. A paint strainer can be used to filter out any debris. Printing with cured debris could potentially cause damage to the tray or distort the printed objects.**
- Deviations from the described manufacturing process may lead to different mechanical and optical properties of the VeriModel™ OS material.

WORKING WITH PRINTED MODELS

- Models can be coated with a separating material such as Whip Mix's Foilcote, Lubritex No. 12, or Gator Die Lube to prevent wax from adhering.
- Models can be duplicated with silicone.
- When models are used to create orthodontic appliances, or denture bases, a non-food grade silicone mold release agent should be used to prevent acrylic from adhering to the surface of the model.
- Models can be used with vacuum form method to create aligners, retainers, and indirect bonding trays.

Danger!

Contains Acrylate, Phosphine Oxides and Methacrylic monomer.

May cause an allergic skin reaction.

Causes serious eye irritation.

Suspected of damaging fertility.

Harmful to aquatic life with long lasting effects.

Precautionary statements:

Obtain special instructions before use.

Do not handle until all safety precautions have been read and understood.

Avoid breathing mists, vapors or spray.

Wash thoroughly after handling.

Contaminated work clothing must not be allowed out of the workplace.

Avoid release to the environment.

Wear protective gloves and eye protection.

IF exposed or concerned: Get medical attention.

IF ON SKIN: Wash with plenty of soap and water.

If skin irritation or rash occurs: Get medical attention.

Take off contaminated clothing and wash it before reuse.

IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do.

Continue rinsing.

Immediately call a POISON CENTER or doctor.

If eye irritation persists: Get medical attention.

Collect spillage.

Store locked up.

Dispose of contents and container in accordance with local and national regulations.

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Whip Mix Corporation ■ 361 Farmington Avenue ■ Louisville, KY USA 40209