Gypsum Bonded Investments

Low-Fusing Crown and Bridge Alloys

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


(Beauty-Cast, Cristobalite and Novocast)

1. Pattern Preparation:

- Attach sprue to the wax pattern using pliable sticky wax.
- Mount patterns on crucible former base.
- Spray with Whip Mix SMOOTHX Wax Pattern Cleaner. Gently blow excess from pattern.

2. Ring Preparation:

- Line casting ring with liner positioned 6–12 mm short on each end.
- Tack down in place with pliable sticky wax.
- Soak lined ring in water for a minimum of one minute, then shake well (5 times).
- Place ring onto former base.

3. Mixing/Investing:

Suggested Water/Powder Ratio:

	50 g	75 g	100 g
Beauty-Cast	15 mL	22.5 mL	30 mL
Cristobalite			
Novocast	17 mL	26 mL	32–36 mL*

* For more expansion, use 32 mL/100 g.
For less expansion, use 36 mL/100 g.

- Rinse bowl out with water — shake out excess.
- Add measured water to bowl then add weighed powder.
- Hand mix for 15 seconds.
- Mechanical mix under vacuum, slow speed (350–450 rpm) for 30–40 seconds.
- Invest immediately under low vibration.

4. Benchset:

- Benchset for a minimum of 30–40 minutes.

5. Burnout:

- Remove crucible former and trim off top of mold.
- Place mold in furnace at room temperature and raise temperature at a rate of 8°C (15°F) per minute to 650°C (1,200°F).
- Hold for 10–15 minutes.

6. Casting:

- Melt alloy with gas/air torch.
- Use only 1–2 turns on centrifugal, broken-arm casting machine to avoid cracking or finning.
- Quench or cool casting according to alloy manufacturer's instructions.

Hygroscopic Technique

Follow Steps 1–3 of Thermal Technique for pattern preparation, ring preparation, mixing and investing.

4. Setting:

- While investment is still fluid, immerse the invested ring under water at 37.8°C (100°F) for 30–40 minutes in HYGROBATH.

5. Burnout:

- Maximum burnout temperature recommended is 515°C (950°F).
- Hold for 30–40 minutes.

Follow Step 6 of Thermal Technique for casting.

Note:

- Use separate mixing bowls and paddles when mixing gypsum and phosphate investments.
- Always use a metal ring when investing gypsum investments.
- Store in a cool, dry place. Gypsum bonded investments will deteriorate if exposed to the atmosphere, especially if the humidity is high.
- For Beauty-Cast Investment, the Hygroscopic Technique provides more expansion than the Thermal Technique.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer). During burnout procedure, investments may give off harmful fumes and adequate ventilation should be used.

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