

slide off if installed correctly. After installing FlexTrack on both fossae, put the upper member back onto the articulator. Your Denar Mark 310 or 320 is now ready to be used as a tracking system.



Figure 3

### **Removal:**

Remove the upper member. There is a small tab on the lower inside part of each FlexTrack. The tab fits into a screw hole on the fossa,


securing it in position. (Figure 3) Pressing with your thumb or finger inside the FlexTrack will release the tab, allowing the FlexTrack to slide off the fossa. Store the FlexTrack in its plastic case when not in use.



RESTORATIVE ORAL HEALTH  
Whip Mix Corporation  
361 Farmington Avenue  
Louisville, KY USA 40209  
Toll-Free: 800-626-5651  
Phone: 502-637-1451  
Fax: 502-634-4512  
www.whipmix.com

© 2012 Whip Mix Corporation  
FN 20001047-F AA 1211

## Denar® FlexTrack Articulator Accessory Instructions



**The FlexTrack provides the flexibility of movement in both open and closed condylar tracking on one articulator. For use on Denar Mark 320 and Mark 310 articulators.**

## The Most Flexible Articulator System: Denar® Mark 320 and Mark 310.

### Closed Tracking Fossae

Delivers consistently accurate movements by keeping the condyles in contact with the fossae at all times during functional movements. Perfect when making guide tables or adjusting eccentric movements since the upper member can not lift away from the lower member, causing inaccurate movements. This will insure the most accurate adjustments of your restorations possible on a semi-adjustable articulator.

### Open Fossae

Traditional “Open Fossae” design allows for diagnostic analysis and planning by permitting separation of condyles and fossae during function. This makes it possible to identify and diagnose occlusal

interferences and other pathologic issues for the patient.

FlexTrack enables you to have both systems in one articulator.

**Note:** FlexTrack does not fit the Denar® Mark 330 articulator.

### Installation:

Remove the upper member by releasing the centric latch and elastic retainers from the shoulder screws. The Denar FlexTrack

consists of a right and left slip on fossa accessory. The opening in the FlexTrack faces forward with the notch toward the fossa. (Figure 1) To install, simply slide the FlexTrack onto the fossa.

You may need to rotate and separate the FlexTrack slightly with your thumb to ease the track until it ‘clicks’ into place. (Figure 2) A small amount of movement in the installed FlexTrack is normal, but the FlexTrack should not easily



Figure 1



Figure 2