VERICAST™ OS

3D PRINT RESIN
INSTRUCTIONS

3D-DRUCKHARZ
ANLEITUNG

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PRODUCT DESCRIPTION

Whip Mix VeriCast™ OS is a 3D print resin for producing dental burnout patterns printed on DLP and LCD technology printers with 385nm and 405nm light sources.

STORAGE

- VeriCast™ OS must be stored in the original packaging at room temperature in a dry area.
- Always keep container tightly sealed and close the container immediately after each use.
- Do not expose VeriCast™ OS to any light source. Minimum amount of light can induce polymerization (solidification).
- Once resin is poured in a tray, it should be stored in a dark area
- **Resin that is stored in a tray must be stirred before each print using a plastic or paper card (similar to a business card), to ensure a homogenous mixture.**

NOTE: Avoid using sharp objects and edges as these can damage the resin tray.

- Do not use the product past the expiration date. The lot number and expiration date are indicated on VeriCast™ OS packaging. In case of questions or concerns, please refer to the lot number.

PROCESSING


- When printing with VeriCast™ OS 3D Print Resin, make sure it is at a temperature between 20°C and 30°C. Extreme high or low temperatures could negatively affect the performance of the resin.
- VeriCast™ OS bottle must be shaken thoroughly for approximately five minutes before use.
- Ensure film on the Build tray is clear of any debris before starting the print.
- Carefully pour VeriCast™ OS into the tray/vat.
- For optimized print settings refer to the latest material file that you can find in the Resource library at <https://whipmix.com/technical-resources/>
- After the build process is finished, post processing is required. If this cannot be immediately completed, leave the printed objects in the printer until ready to post process.

POST PROCESSING

- When the print job is complete, remove the build platform from the printer, then carefully remove the printed objects from the platform.
- Carefully remove all support structures from the printed object.
- **In a well-ventilated area, wash parts for 3 – 5 minutes in alcohol bath.**
- **Printed objects should then be rinsed in a second clean ultrasonic alcohol bath for an additional 5 minutes. Allow parts to dry thoroughly before post curing.**
- To achieve the maximum mechanical strength, post cure in a light curing unit, such as Asiga's Flash from Whip Mix, for 20 minutes. For additional curing units, please refer to the Technical Resources section of the VeriCast OS webpage at www.whipmix.com.
- Look for cured debris in the tray after each print job, especially after a print failure. A paint strainer can be used to filter out any debris. Printing with cured debris could potentially cause damage to the tray or distort the printed objects.
- Deviations from the described manufacturing process may lead to different mechanical and optical properties of the VeriCast™ OS material.

INVESTMENT CASTING / PRESSING

- For optimal results, use Whip Mix ResinVest, which is a phosphate investment made specifically for burning out printed or milled resin patterns. For more information, visit <https://whipmix.com/products/resinvest-investment/>.

Warning!

May cause an allergic skin reaction.
 Suspected of damaging fertility or the unborn child.
 Harmful to aquatic life with long lasting effects.

Precautionary Phrases:

Obtain special instructions before use.
 Do not handle until all safety precautions have been read and understood.
 Avoid breathing mist, vapors or spray.
 Contaminated work clothing must not be allowed out of the workplace.

Avoid release to the environment.
 Wear protective gloves and protective clothing.
 IF exposed or concerned: Get medical advice.
 IF ON SKIN: Wash with plenty of soap and water.
 If skin irritation or rash occurs: Get medical attention.
 Take off contaminated clothing and wash it before reuse.
 Store locked up.
 Dispose of contents and container in accordance with local and national regulations

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Whip Mix Corporation ■ 361 Farmington Avenue ■ Louisville, KY USA 40209