X-20 Chrome Investment

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1. Duplication:

- Duplicate master stone cast with a hydrocolloid.
- Separate stone cast from hydrocolloid mold after 30 minutes and rinse with cold water.
- Lightly blow excess water from inside the mold.

2. Pouring Refractory Cast:

- Pour refractory cast using X-20 investment. For simple metal frames use 50% Special Liquid Concentrate (1 part Liquid, 1 part water). For large metal frames use 100% Special Liquid Concentrate.

11 mL liquid to 100 g powder

- Rinse mixing bowl with water; shake out excess.
- Pour accurately measured liquid into the bowl, then add powder to liquid.
- Hand mix 15 seconds, then mix mechanically at slow speed under vacuum for 30 seconds.
- Vibrate investment into hydrocolloid mold taking care not to trap air bubbles.
- Bench set for 30 minutes.
- Separate refractory cast from duplicating material and air dry for 30 minutes.

3. Trimming Refractory Cast:

- Trim refractory cast with model trimmer removing any sharp angles from cast and taper base of cast to allow outer mold to lock with refractory cast.
- Rinse refractory cast with water to remove any investment particles on surface.

4. Drying Refractory Cast:

- Place refractory cast in cold (room temperature) drying oven.
- Raise to 110°C (230°F) and hold for 1 hour.
- Apply desired surface treatment to refractory cast.
- Drain excess surface treatment if necessary.
- Cast is now ready for wax-up.

5. Investing Wax-Up

- Mount waxed-up cast to base former using Whip Mix PLIABLE STICKY WAX.
- Clean wax-up with Whip Mix SMOOTH EX or other wax pattern cleaner.
- Lightly blow excess cleaner from wax pattern.
- Secure ring around base.
- If using metal ring, line ring with 2" RING LINER leaving liner 6 mm (1/4 inch) short of each end of metal ring.
- Allow lined ring to soak in water for 1 minute, shake ring 5 times to remove excess water, then place ring on base.
- Rinse out mixing bowl with water; shake out excess.
- Pour accurately measured powder into water.

13 mL water to 100 g powder

6. Outer Mold:

- Hand mix 20 seconds to thoroughly incorporate powder into water. Mix mechanically under vacuum for 30 seconds.
- Paint investment into critical areas then vibrate investment into ring taking care not to trap air.
- Bench set for a minimum of 1 hour.
- Trim top of mold and rinse under tap water.

7. Burnout:

- Place into cold oven and raise to 400°C (750°F) using heat rate of 6°C/min (10°F/min).
- Hold for 30–60 minutes, then raise to 900°C (1,650°F) at a heat rate of 10°C/min (20°F/min).
- Hold for 30–60 minutes.

8. Casting:

- Cast, then let mold cool to room temperature.
- Devest and sand blast in usual manner.

Note:

- Recommended Liquid: Special Liquid Concentrate
- For changes in expansion rates please call Whip Mix Technical Support.
- Always use separate mixing bowls for phosphate and gypsum investments.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

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