# Prestobalite

Gypsum Investment For Rapid Casting of  
Low-Fusing Crown and Bridge Alloys

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
## Technique

### 1. Pattern Preparation:

- Attach sprue to pattern in usual manner and mount on crucible former.
- Clean pattern with a surface tension-reducing agent such as Whip Mix SMOOTHTEX.
- Dry thoroughly with a gentle blast of air.

### 2. Ring Preparation:

- Use PLIABLE STICKY WAX to hold 1 layer of Whip Mix RING LINER in place.
- Leave liner short on both ends of the ring by 6 mm.
- Soak in water for a minimum of one minute to saturate liner then shake off excess water from ring.
- Do not blot or touch wet liner.

**Note:** Prestobalite must be used with a metal casting ring. We also recommend a paper ring liner as evidence shows that ceramic liners contain fibers which may be hazardous to your health. If a ceramic liner is used, we recommend soaking the liner in Whip Mix Smoothtex debubblizer.

### 3. Water/Powder Ratio: 19 ml/60 g

- Expansion may be controlled by varying the water/powder ratio.  
For more expansion, use 18 ml/60 g  
For less expansion, use 20 ml/60 g

### 4. Mixing and Investing:

- Wet bowl and shake out excess.
- Pour measured water into bowl and add investment powder.
- Stir and mix with spatula until powder is fully incorporated.
- Mechanical mix under vacuum, slow speed (350-450 RPM) for 60 seconds.
- Invest ring.

### 5. Setting:

- Benchset for at least 30 minutes.
- Trim top of mold.
- Remove crucible former.

### 6. Burnout:

- Place casting ring in burnout furnace preheated to 650°C (1200°F) for a minimum of 30 minutes. Longer times may be needed with multiple rings.
- When using plastic or vinyl sprues or posts, it is best to start with a cold oven and raise the temperature to 315°C (600°F) and hold for 30 minutes. Then raise temperature to 650°C (1200°F).

### 7. Casting:

- When casting with a centrifugal casting machine, 2 turns is sufficient. Too much force when casting may cause "finning or flash" on the metal casting.

### 8. Recovery of Casting:

- Let the casting cool for 2 minutes.
- Then quench in water and devest.

### Warning

Investments contain free silica – DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer). During burnout procedure, investments may give off harmful fumes and adequate ventilation should be used.

### Helpful Hints

**For larger castings, please note the following suggestions:**

- Use less water
- Increase mixing time
- Warm the water before mixing

**For smaller castings, please note the following suggestions:**

- Use more water
- Decrease mixing time
- Cool the water before mixing

### Note:

- This investment is for alloys with a melting point below 1080°C (1976°F).
- Gypsum bonded investments are subject to deterioration when exposed to the atmosphere, especially in high humidity climates.

---

#### Physical Properties:

Water/Powder Ratio	32 mL/100 g
Working Time	5–7 minutes
Setting Time (ADA Method)	15 minutes
Ready for Burnout	30 minutes
Setting Expansion	0.6%
Thermal Expansion – 700°C	1.1%
Compressive Strength, Wet (2 hrs.)	500 psi (3.4 MPa) minimum

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