



more reasons. one source

Whip Mix Corporation • 361 Farmington Avenue • P.O. Box 17183 • Louisville, Kentucky 40217 USA • phone: (502)637-1451 • fax(502)634-4512 • www.whipmix.com

## Whip Mix Introduces the MainStay™ II

**Louisville, KY: November 2009** – Whip Mix adds a dual dowel pin system, the MainStay™ II, to the already popular MainStay™ line of products. The MainStay™ II is a uniquely designed system for fabricating crown and bridge working models with sectioned, removable dies. Its protective sleeve and tapered, dual pin design prevents reverse placement of sectioned dies. In addition the protective sleeve grips the dowel pin creating positive indexing while protecting the base stone material from abrasion. The pins are compatible with most pin setting machines, therefore requiring no change in procedure. MainStay™ II's rubber end cap stabilizes quadrant models prior to base stone application and makes it easy to expose the dowel pin without the use of a dowel pin bur. MainStay™ II is available as a box of 1000 pins, sleeves and rubber end caps. For a limited time, Whip Mix is offering a special introductory promotion, Buy 3 boxes, get 1 FREE.

For more information contact Whip Mix Corporation at 800-626-5651 or [click here](#) for a FREE sample.

*Founded in 1919, Whip Mix manufactures Hanau™, Denar® and Whip Mix articulators, dental gypsums, investments, waxes, furnaces, ovens, mixing equipment and model trimmers for use in dental labs and dental offices. Whip Mix Dental Technology Solutions offers the Vericore® System, which utilizes a state-of-the-art CAD/CAM technology to mill custom Zirconia copings bridges and abutments and Titanium implant, supported milled bars. Whip Mix, headquartered in Louisville, Ky, is ISO9001-2000 certified and sells its products through distributors in over 80 countries and also maintains an office and distribution center in Dortmund, Germany. For more information visit [www.whipmix.com](http://www.whipmix.com).*

