



more reasons. one source

Whip Mix Corporation • 361 Farmington Avenue • P.O. Box 17183 • Louisville, Kentucky 40217 USA • phone: (502)637-1451 • fax(502)634-4512 • www.whipmix.com

WHIP MIX INTRODUCES THE Denar[®] MARK 300 SERIES ARTICULATORS

Louisville, Kentucky: March 2009 – Building on a foundation known for its comprehensive systems and strong engineering, Whip Mix is proud to introduce the next generation of articulators, the Denar[®] Mark 300 Series Articulators. The Mark 300 Series offers interchangeability among the instruments and are factory set to within 20 microns of accuracy.

The Mark 300 Series is comprised of the Mark 330, Mark 320 and Mark 310 Articulators. The Mark 330 is semi-adjustable with a condylar inclination of 0-60 degrees, adjustable immediate side shift of 0-2mm and adjustable side shift of 0-15 degrees. The Mark 320 is also semi-adjustable with an adjustable condylar inclination of 0-60 degrees and progressive side shift fixed at 15 degrees. The Mark 310 has fixed settings with a condylar inclination of 25 degrees and a progressive side shift of 15 degrees. The Mark 300 Series has a new magnetic mounting system, which offers 3 points of contact to ensure stability. In addition, the Mark 300 Series is compatible with the Slidematic Facebow. There is also a positive centric latch which allows the upper and lower members to be separated or positively locked together.

The new series enhances the field of articulators as we know them – they're easier to use, have a new sleek design, helps you obtain predictable results and are factory set to unsurpassed accuracy. For more information on the Mark 300 Series Articulators please visit our website at www.whipmix.com or call 1-800-626-5651.

Founded in 1919, Whip Mix manufactures Hanau[™], Denar[®] and Whip Mix articulators, dental gypsums, investments, waxes, furnaces, ovens, mixing equipment and model trimmers for use in dental labs and dental offices. Whip Mix Laboratory Services Division offers the Vericore[™] System which utilizes a state-of-the-art CAD/CAM technology to mill custom zirconia copings and bridges. Whip Mix is ISO9001-2000 certified and sells its products through distributors in over 80 countries and maintains an office and distribution center in Dortmund, Germany. For more information visit www.whipmix.com.



Mark 330



Mark 320



Mark 310