





ON-GOING WHIP MIX PRINT RESIN/ 3D PRINTER QUALIFICATION

| QUALIFICATION COMPLETE = ▲ ONGOING QUALIFICATION PROCESS = ● | VERIMODEL OS* | | | | VERITRAY* (385/405nm) | VERICAST* OS (385/405nm) | VERIGUIDE* OS CLEAR (385/405nm) | VERISPLINT** OS CLEAR 510K Cleared (385nm) | DENTCA*** DENTURE/ TRY-IN/ TEETH | DENTCA*** CROWNS & BRIDGE | SURGICAL GUIDE FOR CARBON | KEYSPLINT*** SOFT RESIN | BEGO VARSEOSMILE CROWN PLUS AND BEGO VARSEOSMILE TEMP |
|--|----------------------|----------------------|---------------------|--------------------------------|--------------------------|--------------------------------|---------------------------------------|---|---|---------------------------------|---------------------------------|-------------------------------|--|
| | IVORY (385/405nm) | WHITE (385/405nm) | GREY (385/405nm) | GOLDEN BROWN (385/405nm) | | | | | | | | | |
|  WHIP MIX VERIEKO LCD (405nm) | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | | ▲ | ▲ |
|  WHIP MIX VERIBUILD LCD (405nm) | ▲ | ▲ | ▲ | ▲ | | ▲ | ▲ | | ▲ | ▲ | | | ▲ |
|  ASIGA MAX (385nm/405nm) | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | | ▲ | ▲ |
|  ASIGA PRO 2 75 (385nm/405nm) | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | | | |
|  ASIGA PRO 4K (385nm/405nm) | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | ▲ | | ▲ | ● |
|  ACKURETTA FREESHAPE 120 (405nm) | ▲ | ▲ | ▲ | ▲ | | ▲ | ▲ | | | | | ▲ | ▲ |
|  ACKURETTA DIPLO (405nm) | | ▲ | ● | ▲ | | ▲ | ▲ | | | | | | |
|  ANYCUBIC (405nm) | | ▲ | | ▲ | | ▲ | ▲ | | | | | | |
|  FLASHFORGE HUNTER (405nm) | | ▲ | | ● | ▲ | ▲ | ▲ | | | | | | |
|  MICROLAY VERSUS (385NM/405NM) | ● | ▲ | | ▲ | ▲ | ▲ | ▲ | ● | | | | | |
|  NEW PRO 3D NP1 (385nm/405nm) | ● | ▲ | | ▲ | ▲ | ▲ | ● | ● | | | | | |
|  CARBON 3D M1, M2, M2D | | | | | | | | ● | | | ▲ | | |
|  MAKEX M-DENTAL | | ▲ | | | | ▲ | ▲ | ● | | | | | |
|  MAKEX M-ONE PRO | | | | ▲ | ▲ | | | | | | | | |
|  PHOTOCENTRIC LC DENTAL | | | | ▲ | ▲ | ▲ | ▲ | ● | | | | | |

* Available globally • ** Available in the U.S. only • *** Available in the U.S. and Canada only
Don't see your printer? Contact us at www.whipmix.com/contact-us Visit www.whipmix.com for updates.