

VERILUX CURING UNIT

MAKING LIGHT WORK OF POST-PROCESSING.



THE WHIP MIX VERILUX IS HERE!

Use the VeriLUX light curing unit as a very affordable option for the effective polymerization of light-cured resins.

Meeting the needs of both dentists and labs, the VeriLUX curing unit is small, yet provides a sufficient-size curing chamber. Best of all, it is fast and super-easy to use.


*Whip Mix VeriLUX
Item #79435*

- Easy to use — Preloaded validated curing time database for wide range of resins with the turn of a dial.
- Accuracy and consistency — Preloaded curing settings eliminate entering the wrong time or intensity.
- Superior Performance — Smart UV-light source management system for maximized energy efficiency.
- High efficiency — Cures most resins within a short time and, in most cases, there's no need to cure both sides.
- Versatility — Add personalized settings by changing the time, power, frequency.
- Very affordable!

502.637.1451 • whipmix.com



WHIP MIX VERILUX SPECIFICATIONS



Curing Time	1 second to 20 minutes
Light Intensity	190 mW/cm ²
Lamp Wavelength	UVA + UVV (365-405nm)
Product Dimensions	200 x 220 x 160 mm
Curing Chamber Dimensions	140 X 165 X 90 mm
Power Supply	100-240 V AC, 50-60 Hz
LED Power Output	75 W
Fuse	2A
Net Weight	3Kg
Connectivity	SD Card

YOUR VERILUX CURING UNIT COMES WITH:

- SD card reader for updating and saving material settings and curing times.
- One-year warranty
- Glass plate
- Quick Start Guide
- Power Cable

"We love testing out new products for Whip Mix! The VeriLUX Curing Unit was no exception. Honestly, I was shocked at how tiny, but mighty this curing unit is. We have a couple of different types of units and I can honestly say, I was so pleased with its efficiency and power. The quickness of the cure makes it a true value to same day dentistry. I'm pleased to add this unit to our army of Whip Mix products that we love!"

Brandon Stapleton, DMD, MSD, iDENTity Dental Studio — Lexington, KY

VISIT WHIPMIX.COM OR CALL 502.637.1451 TO LEARN MORE ABOUT HOW VERILUX CAN STREAMLINE YOUR PRINTING WORKFLOW.



WHIP MIX CORPORATION
361 Farmington Avenue • Louisville, KY USA 40209
P: 502-637-1451 • whipmix.com