


# DENAR®/HANAU™ PARTS LIST

## DENAR/HANAU VERIFICATION GAGE

20011650

PART #: 20011681  
EFFECTIVE: 3-1-2010



Key #	Part #	Description
1		Upper Assembly*
2		Lower Assembly*
3	20011656	Size Pin, Dia. .246
4	20011655	1/4" x 1" Dowel Pin

Key #	Part #	Description
5	3037	Screw, Socket Head Cap
6	3178	Knob
7	3138	Screw, Socket Head Cap**

\* Factory Repair Required  
\*\* With #6 used for Upper Member of Wide-Vue



RESTORATIVE ORAL HEALTH  
Whip Mix Corporation  
361 Farmington Avenue  
Louisville, KY USA 40209  
Toll-Free: 800-626-5651  
Fax: 502-634-4512  
www.whipmix.com