


Occlusal Rims – Without The Flame!



Electric Occlusal Rim Former

Whip Mix's Electric Occlusal Rim Former offers the technician, the dentist or the dental student the ability to create the ideal shape and size of occlusal (bite) rims made out of base plate wax **without a flame**. With a variable temperature knob, the small unit cleanly melts any wax. In addition, the paddle holder also comes complete with a cover to prevent burn injuries.

- Electric Heating
- Paddle temperature: 0°F – 220°F / 0°C – 104.4°C
- Dimensions: 75 mm W x 230 mm D x 105 mm H
- Power Usage: 100 W
- 2-Year Warranty



Whip Mix Corporation
800-626-5651 502-637-1451 www.whipmix.com

Connect With Us!

