



USING VERICAST OS 3D PRINT RESIN MEANS ACCURATE CASTINGS & PRESSED CERAMICS.

- Prints detailed crowns, bridges, substructures, and RPD frameworks accurately, smoothly, is durable, and burns out cleanly.
- Red color makes it easy to see detail.
- Ideal for investment casting and ceramic pressing.
- Works with LCD and DLP printers in both 385nm and 405nm wavelengths.



For optimal results, use ResinVest phosphate investment made especially for burning out printed or milled resin patterns.



502.637.1451 • 800.626.5651 • whipmix.com

DENTAL TECHNOLOGY SOLUTIONS