



Flask Band

Older and mismatched flasks often have spaces between the flask parts



Flowable stones leak from these sites, creating a mess



An incomplete flask pour



Flask Band



The silicone band and block seals the openings between the upper and lower flask parts



Resulting pour is smooth & accurate with a minimum amount of flash



- Flask Band reduces leakage of flasking stone in older and mismatched flasks
- Neater/Cleaner flasking procedure
- Extends the life of your brass flasks by improving the seal of the upper and lower connecting surface
- Enables the use of flowable Flasking Stone resulting in a smoother and more accurate reproduction

Item # 33977 - 3 Flask Bands

To order call: 800.626.5651