

NEW!

Flexible Vac-U-Mixer

FLEXIBLE BOWL FOR MIXING INVESTMENTS AND GYPSUMS UNDER VACUUM

Few things frustrate model technicians more than stone that has hardened in an expensive bowl. They must scrape the bowl manually, often damaging it beyond use. Many times, the hardened stone can't be removed at all.

Enter the revolutionary Flexible Vac-U-Mixer for the Vacuum Power Mixer, Combination Unit and VPM2. With just a firm squeeze, the technician can crack and release a majority of the hardened stone without damaging the bowl.

THE FLEXIBLE VAC-U-MIXER ADVANTAGE

- Easy clean up
- First ever flexible bowl for vacuum mixing – Higher impact resistance than a standard mixing bowl
- Visual indicators of the max fill line – Prevents overfilling the bowl which can result in damage to mixing unit
- Larger opening for easier hand access – Contoured for easy spatula access
- 350 g bowl – Lid & Paddle Assembly included

Patent Pending

FIRM SQUEEZE
RELEASES HARDENED STONE



VACUUM MIXING



EASY ACCESS



Item#	Description
15665	#6720 Flexible Vac-U-Mixer for use with the Vacuum Power Mixer and Combination Unit (Lid & Paddle Assembly included)
29352	#6785 Flexible Vac-U-Mixer for use with the VPM2 (Lid & Paddle Assembly included)

For more information, call 800-626-5651 x2 or visit our website at www.whipmix.com.



more reasons. one source