


WHIP MIX DTS PRINTING SERVICES

**Our Dental Technology Solutions outsource services now include
3D printing of models, surgical guides, and castable patterns.**



HOW TO WORK WITH US

Dental Technology Solutions accepts all .stl files. If you have a 3Shape scanner, we can add the Whip Mix material files to your 3Shape selections. If you do not have a 3Shape scanner, upload your file to our web site.

For more information contact Customer Service at 800-626-5651, option 2.

whipmix.com



DENTAL TECHNOLOGY SOLUTIONS

WHIP MIX 3D PRINT SERVICES

NEED 3D PRINTED PARTS? TOO MUCH PRODUCTION TO HANDLE? 3D PRINTER TEMPORARILY DOWN? OVERWHELMED? SEND US YOUR FILES — SIMPLY VISIT OUR WEBSITE.

Whip Mix is your perfect partner when you need printed models, surgical guides or even printed C&B patterns for casting or ceramic pressing with a 96-hour turnaround time. We're easy to work with and very affordable.



PRINTED MODELS

Accurate, smooth and detailed Vericore Printed models are available in Full arch, Quadrant, or even POM (Peace of Mind) configurations, with or without removable dies. Printed in ivory color.

VeriModel OS resins produce 3D printed models with reliable accuracy, fine detail and clean, smooth surfaces.



SURGICAL GUIDES

The Surgical Guides are made from VeriGuide OS Print Resin, a Class I, CE certified resin material developed for the production of high precision surgical guides for use in dental implant surgery.

VeriGuide OS™ has been successfully tested for biocompatibility and meets all mechanical properties requirements. The resin produces highly accurate surgical guides with a smooth, high quality surface and provides an excellent fit for sleeve insertion. The material can also be sterilized by autoclave, steam, ethylene oxide, or gamma rays without affecting dimensional stability, physical properties and biocompatibility.



RPD/C&B PATTERNS FOR CASTING/CERAMIC PRESSING

VeriCast OS adds a Burnout Resin.

- Prints detailed crowns, bridges, substructures, and RPD frameworks accurately, smoothly, is durable, and burns out cleanly.
- Red color makes it easy to see detail.
- Ideal for investment casting and ceramic pressing.
- Works with LCD and DLP printers in both 385nm and 405nm wavelengths.

3D PRINT SERVICES

79085	Vericore Printed Model — Full Arch OS Ivory
79086	Vericore Printed Model — Quadrant OS Ivory
79087	Vericore Printed Model — POM Segment OS Ivory

ADDITIONAL SERVICES

79500	Printed Model Service Fee for Model Builder for all scans
-------	---

SURGICAL GUIDES (Translucent)

79075	Surgical Guide
-------	----------------

CASTABLE PATTERNS (For C & B and partial frameworks)

79060	Castable Pattern — Crown (Up to 6 units)
79076	Castable Pattern — RPD, Long Span Bridge



WHIP MIX CORPORATION
 361 Farmington Avenue • Louisville, KY USA 40209
 P: 502-637-1451 • TF: 800-626-5651 • W: whipmix.com