

# New! Denar<sup>®</sup> Best-Bite<sup>™</sup> Discluder



*Help patients feel their best*




*It's hard to look good when you're in pain. In a landmark study, Dr. Jannet Travell, MD found that the pain caused by conflict between the teeth and jaw joints can radiate and refer to other areas of the head, neck and face with varying intensity and frequency. These symptoms of occlusal interference are commonly confused with migraine, tension, cluster and stress headaches.*

*Call*

**800-626-5651**

today to find out how easy the Denar<sup>®</sup> Best-Bite<sup>™</sup> Discluder is to use.

Every day at least one patient visits your dental office with head, neck and facial pain. Now you can diagnose whether occlusal interference is the source of patients' pain quickly, easily and reliably with the Denar<sup>®</sup> Best-Bite<sup>™</sup> Discluder.



### ***Easy to Use***

Simple 4 step diagnostic process.

### ***Rapid Results***

There will either be noticeable pain relief or not within a few minutes.

### ***Safe and Effective***

When used as instructed, this is a safe and quick tool to identify whether occlusal interferences are causing muscle soreness.

## BEST-BITE<sup>™</sup>



The Denar<sup>®</sup> Best-Bite<sup>™</sup> Discluder is an indispensable part of any dental practice.



RESTORATIVE ORAL HEALTH

Whip Mix Corporation  
361 Farmington Avenue  
Louisville, KY USA 40209  
Toll-Free: 800-626-5651  
[www.whipmix.com](http://www.whipmix.com)

# BEST-BITE™



***Determining whether pain is from muscle spasms due to occlusal interference is a simple four step process:***

**1**

Check for clenching and grinding of teeth.

**2**

Check for occlusal conflict between teeth and jaw joints.

**3**

Check for corresponding muscle soreness.

**4**

Use the Best-Bite Discluder to see if the pain from sore muscles is relieved.



***Available in two convenient kits:***

**8 Patient Introductory Kit (Item #20000011)** – Treats up to 8 patients and includes the complete Best-Bite educational package including educational videos on CD, detailed doctor and patient instructions, safety straps, and storage cases.


**24 Patient Complete Refill Kit (Item #20000012)** – This kit is recommended for those who have already purchased the Introductory Kit. It treats up to 24 patients and includes detailed doctor and patient instructions, safety straps, and storage cases.



RESTORATIVE ORAL HEALTH  
Whip Mix Corporation  
361 Farmington Avenue  
Louisville, KY USA 40209  
Toll-Free: 800-626-5651  
[www.whipmix.com](http://www.whipmix.com)

***What does the Best-Bite™ Discluder do***

**?**



*Centered jaw and teeth fit.*



*Centered jaw and teeth do not fit.*



*Teeth pull jaw off center leading to muscular pain.*



*Best-Bite™ Discluder aligns the jaw without interference from the teeth, allowing the jaw joint to return to its natural position.*