

# VERIEKO PRINTER

PRICE. SPEED. PERFORMANCE.  
NO COMPROMISE.



## THE WHIP MIX VERIEKO IS HERE!

This outstanding compact printer delivers a printing speed of up to 4.5 cm/hour with exceptional accuracy, due to its 95+% LED light uniformity.

VeriEKO is ready to prove its printing power right out of the box!

*Whip Mix VeriEKO  
Item #79401*

- 
- 54 calibrated LEDs means curing consistency
  - Ultra-accurate, smooth prints with 49µm pixel-size
  - Anti-drip LCD panel protection
  - Easy-to-use touchscreen
  - Fail-safe VAT screen
  - Minimal learning curve
  - Easy replacement parts
  - Industry-leading 2-year warranty
  - Intuitive and dental-focused 3D slicing software

**502.637.1451 • [whipmix.com](http://whipmix.com)**



DENTAL TECHNOLOGY SOLUTIONS

# WHIP MIX VERIEKO SPECIFICATIONS



<b>Technology</b>	Monochrome LCD
<b>Printing Size</b>	128 x 80 mm
<b>Pixel Size</b>	49 micron
<b>Printing Speed</b>	Up to 4.5 cm/hr
<b>LCD</b>	Up to 10,000 operation hours
<b>Backlight</b>	54 LEDs

## COMPATIBILITY AND PERFORMANCE

The Whip Mix VeriEKO printer and specialty resins are ideal for a long list of dental applications such as permanent crown and temporary crowns and bridges, veneers, inlay and onlay, models with removable dies, splints, and surgical guides. We are continually working to qualify more resins.

## YOUR VERIEKO PRINTER COMES WITH:

- 2 Extra Vats
- 1 Wi-Fi dongle
- 1 Extra LCD Print Screen
- 1 Hex key
- 1 Ethernet cable
- 1 Power cable
- FREE Alpha 3D Software License

*"Overall I am quite impressed with this little printer...it's outstanding. Very simple to use, very clean, precise, and smooth. It's beautiful visually with a very clean design inside."*

*Michael Scherer, DMD, MS — Sonora, CA*

**VISIT WHIPMIX.COM OR CALL 502.637.1451 TO LEARN MORE ABOUT HOW VERIEKO CAN STREAMLINE YOUR PRINTING WORKFLOW.**



**WHIP MIX CORPORATION**  
361 Farmington Avenue • Louisville, KY USA 40209  
P: 502-637-1451 • [whipmix.com](http://whipmix.com)