

# VERIPULSE IOS

PRECISION YOU EXPECTED.  
SIMPLICITY YOU DID NOT.



**Whip Mix VeriPULSE**  
**Item # 70675**

The VeriPulse IOS is a solution that is ready to use right out of the box that offers premium user experience, that requires little training, implements easily in the daily workflow and low learning curve which makes it the perfect choice for daily use by dental professionals worldwide.

With next-generation software utilizing new and innovative technology, the VeriPulse IOS greatly increases scan speed and accuracy while adding even more features to support lab, restorative, and implant workflows.

- Acquisition PC — Veripulse IOS comes complete with an optimized PC and pre-installed software.
- Truly Ergonomic — Compact design, weighs only 150 grams, autoclavable 360° rotating anti-fogging tip for optimum scanning angle
- Truly Open System — Compatible with most dental CAD systems providing flexibility in the lab's choice of CAD system. Output format STL, PLY and OBJ.
- High Efficiency — Fast scanning, full arch scan time < 2 minutes
- Accuracy — 22µm quadrant accuracy, HD image capture and color scanning, with no powder or spray required
- Cloud-based Sharing — Patient information, including scans and all other files, is stored and shared with your lab using our cloud-based sharing platform, WhipMixCloud
- Color Scanning 24-bit (8-bit per channel)
- VeriPULSE features a 2-year warranty
- No Fees. Ever.



**502.637.1451 • [whipmix.com](http://whipmix.com)**

DENTAL TECHNOLOGY SOLUTIONS

# VERIPULSE IOS SOFTWARE

## WhipMix Clinic -> WhipMix Cloud



### **Case Creation**

Search for or create a patient case and manage the order in just a few steps.



### **Scanning**

A fast, intuitive process performed by the dentist or the dental assistant.



### **Case Sharing**

Sharing the digital impression with your lab of choice.



### **Lab Integration**

The software is an open system, compatible with all CAD systems.

## **WHIPMIX CLINIC**

WhipMix Clinic Software is a highly intuitive scan and order management software package that is included with the VeriPulse IOS. In addition to order management and scanning, the system allows you to share information with the lab via a cloud file-sharing solution. The software is highly intuitive, with a fast learning curve that provides fast and easy collaboration with the lab. WhipMix Cloud will provide the lab you have configured a connection to, and to which you have sent the order from HeronClinic, with secure access to all patient case files associated with the order.

Using WhipMix Clinic ensures:

- A user-friendly, easy-to-navigate interface
- A clear workflow with a low learning curve for staff
- You can meet all of your restorative, orthodontic, and implant needs
- You can easily create or browse patient cases, scan, and share with your lab

**VISIT [WHIPMIX.COM](http://WHIPMIX.COM) OR CALL 502.637.1451 TO LEARN MORE ABOUT HOW VERIPULSE IOS CAN STREAMLINE YOUR OFFICE WORKFLOW.**



**WHIP MIX CORPORATION**

**361 Farmington Avenue • Louisville, KY USA 40209**

**P: 502-637-1451 • [whipmix.com](http://whipmix.com)**