

3Shape TRIOS MOVE+

Frequently Asked Questions

For reseller use only



Product description

What is the new TRIOS® MOVE+?

TRIOS MOVE+ is a new larger TRIOS monitor. The workspace is 36% bigger than TRIOS MOVE - 15.6" vs. 13.3"- and it comes with a USB port at the back of the screen.

What is the price of TRIOS MOVE+?

Please see TRIOS price list.

Is TRIOS MOVE+ available for all TRIOS scanner types?

Yes.

Is MOVE+ heavier than the original TRIOS MOVE?

No. To provide great ergonomics, MOVE+ has been designed to be the same weight 14 kg.

Which graphics card does TRIOS MOVE+ come with?

MOVE+ has the same graphics card as TRIOS MOVE.

Which type of USB port does TRIOS MOVE+ have?

The USB port on the rear side of the screen is a USB 2.0.

Can MOVE+ be connected to an external screen?

Yes, the MOVE+ supports video output through the USB port on rear side of monitor using the JUA350 USB to HDMI from j5create.

Does TRIOS MOVE+ come with a Stylus Pen?

Yes, a new high-quality Stylus Pen is included in the package.

Is it possible to upgrade TRIOS MOVE to TRIOS MOVE+ by changing the screen only?

No.

Is there a trade-up option for TRIOS MOVE to TRIOS MOVE+?

No.

How is the dongle and subscription handling when TRIOS MOVE+ is ordered as standalone?

Please include the dongle number of the TRIOS software subscription when ordering MOVE/MOVE+ standalone to maintain a seamless subscription.

What are the technical specifications for MOVE+?

More technical details can be found in the Info Pack on Partner Portal.

Availability

Where is MOVE+ available?

MOVE+ is available in all countries where the product has been registered. Please see the TRIOS Registration Status on the Partner Portal.