



COMPATIBLE PRINTERS WITH WHIP MIX VeriRESIN 3D PRINT MATERIAL

Validation Complete Validation in Progress	VeriMODEL All Purpose*			VeriMODEL Ortho*	VeriTRAY**	VeriCAST*	VeriGUIDE**	VeriSPLINT***	VeriGUM*	VeriDENT Try-In**
	IVORY	GOLDEN BROWN	GREY	WHITE	TURQUOISE	RED	CLEAR	CLEAR	PINK	A2
VERIEKO (405nm)	✓	✓	✓	✓	✓	✓	✓		✓	✓
VERIBUILD (405nm)	✓	✓	✓	✓		✓	✓			✓
ASIGA MAX (385/ 405nm)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 2 75 (385nm)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ASIGA PRO 4K (385nm)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
ACKURETTA SOL (405 nm)	✓	✓	✓	✓	✓	✓	✓		✓	✓
ACKURETTA FREESHape 120 (405nm)	✓	✓	✓	✓		✓	✓			
ACKURETTA DIPLO (405nm)		✓	▶▶	✓		✓	✓			
RAPIDSHAPE D20+ / D30+	▶▶	✓	✓		✓		✓		✓	✓
RAPIDSHAPE D50+	✓	✓	✓		✓				✓	✓
STRUMANN P20+ /P30+	▶▶	✓	✓		✓		✓		✓	✓
FLASHFORGE HUNTER (405nm)		✓		▶▶		✓	✓			
MICROLAY VERSUS (385/405nm)	▶▶	✓		✓		✓	✓	▶▶		
NEW PRO 3D NP1(385/405nm)	▶▶	✓		✓		✓	▶▶	▶▶		
CARBON 3D M1, M2, M2D							✓			
SHINING3D AccuFab-L4D Dental (405nm)	✓	✓	✓		✓	✓	✓		✓	✓
PHROZEN SONIC 4K (405nm)					✓				✓	
PHROZEN SONIC MINI (405nm)					✓				✓	
PHROZEN SONIC XL 4K(405nm)	✓	✓	✓							
RAYSHAPE SHAPE 1+DENTAL (405nm)	✓	✓	✓	✓	✓	▶▶	✓		✓	✓
RAYSHAPE EDGE E1 (405nm)	✓	✓	✓	✓	✓	▶▶	✓		✓	✓
UNIZ NBEE (405nm)	✓	✓	✓	✓	▶▶	▶▶	▶▶		✓	▶▶
MAKEX M-DENTAL		✓				✓	✓			
MAKEX M-ONE PRO				✓						

*Available globally / ** Available in U.S., Canada and Latin America only / *** Available in U.S. and Canada only